

BD—Your partner in immunization


2020–2021 product recommendations

Buy six, get one free



BD Eclipse™ Safety Needles and Syringes

The safety mechanism activates with a single finger stroke immediately after injection—no hard surface is required. The wide, textured finger pad provides safe activation for the thumb or forefinger.



BD Integra™ Safety Needles and Syringes

The needle retracts into the syringe, which eliminates exposure to the needle after activation. The clinician immediately activates the retracting mechanism after the needle is withdrawn from the patient.



BD SafetyGlide™ Safety Needles and Syringes

BD SafetyGlide™ products feature single-handed activation that allows for ease of use and minimal training.

BD Activation-Assist™ technology minimizes the finger motion required for activation and keeps fingers behind the needle tip.



BD—Canada, Mississauga, ON, L5N 0B3
905.288.6000

bd.com

BD, the BD Logo, Eclipse, Integra, Luer-Lok and SafetyGlide are trademarks of Becton, Dickinson and Company or its affiliates.
© 2020 BD. All rights reserved.



BD 2020–2021 product guidelines for vaccination administration

BD Safety Needles and Syringes

Configuration	BD product code	BD product description	Packaging
BD Eclipse™ Safety Needles and Syringes			
Needle only	305759	25 G x 5/8 in. BD Eclipse™ Needle	100/shelf pack, 1,200/case
Needle only	305761	25 G x 1 in. BD Eclipse™ Needle	100/shelf pack, 1,200/case
Needle only	305762	23 G x 1 in. BD Eclipse™ Needle	100/shelf pack, 1,200/case
Needle only	305763	22 G x 1.5 in. BD Eclipse™ Needle	100/shelf pack, 1,200/case
Needle only	305767	25 G x 1 1/2 in. BD Eclipse™ Needle	100/shelf pack, 1,200/case
Combination	305780	1 mL BD Luer-Lok™ Syringe with BD Eclipse™ Needle 25 G x 5/8 in.	50/shelf pack, 300/case
Combination	305781	3 mL BD Luer-Lok™ Syringe with BD Eclipse™ Needle 25 G x 5/8 in.	50/shelf pack, 300/case
Combination	305787	3 mL BD Luer-Lok™ Syringe with BD Eclipse™ Needle 25 G x 1 in.	50/shelf pack, 300/case
Combination	305783	3 mL BD Luer-Lok™ Syringe with BD Eclipse™ Needle 22 G x 1.5 in.	50/shelf pack, 300/case
BD SafetyGlide™ Safety Needles and Syringes			
Needle only	305901	25 G x 5/8 in. BD SafetyGlide™ Needle	50/shelf pack, 500/case
Needle only	305916	25 G x 1 in. BD SafetyGlide™ Needle	50/shelf pack, 500/case
Needle only	305902	23 G x 1 in. BD SafetyGlide™ Needle	50/shelf pack, 500/case
Needle only	305900	22 G x 1 1/2 in. BD SafetyGlide™ Needle	50/shelf pack, 500/case
Combination	305903	1 mL BD Luer-Lok™ Syringe with BD SafetyGlide™ Needle 25 G x 5/8 in.	50/shelf pack, 400/case
Combination	305904	3 mL BD Luer-Lok™ Syringe with BD SafetyGlide™ Needle 25 G x 5/8 in.	50/shelf pack, 400/case
Combination	305924	3 mL BD Luer-Lok™ Syringe with BD SafetyGlide™ Needle 25 G x 1 in.	50/shelf pack, 400/case
Combination	305906	3 mL BD Luer-Lok™ Syringe with BD SafetyGlide™ Needle 22 G x 1.5 in.	50/shelf pack, 400/case



Configuration	BD product code	BD product description	Packaging
BD Integra™ Safety Needles and Syringes			
Combination	305269	3 mL BD Integra™ Syringe with BD Integra™ Needle 25 G x 5/8 in.	100/shelf pack, 400/case
Combination	305270	3 mL BD Integra™ Syringe with BD Integra™ Needle 25 G x 1 in.	100/shelf pack, 400/case
Combination	305272	3 mL BD Integra™ Syringe with BD Integra™ Needle 22G x 1.5 in.	100/shelf pack, 400/case
BD Syringe			
Syringe only	309628	1 mL BD Luer-Lok™ Syringe	100/shelf pack, 800/case
Syringe only	309657	3 mL BD Luer-Lok™ Syringe	200/shelf pack, 800/case

Promotion is limited to item numbers listed as purchased by nonacute care customers. Qualifying product purchases must be made between August 1, 2020–December 31, 2020. Invoices must show purchase of six or more items from one of the BD catalogue numbers listed above to receive one free additional product of the same catalogue number. Free offers are limited to 10 free boxes (*needles and syringes*) per customer location. May not be used with any other discount, coupon or offer. BD reserves the right to rescind, revoke or amend this offer without notice.

BD–Canada, Mississauga, ON, L5N 0B3
905.288.6000

bd.com

BD, the BD Logo, Eclipse, Integra, Luer-Lok and SafetyGlide are trademarks of Becton, Dickinson and Company or its affiliates.
© 2020 BD. All rights reserved. (4928/0620)





Be prepared and save

Act now

Take advantage of a limited-time offer from BD. Get **FREE** product when you purchase BD needles, syringes and sharps collectors.

Buy your products all at once or any time between August 1, 2020–December 31, 2020.

Buy 6, get 1 free

Rules and eligibility:

1. Promotion runs between August 1, 2020–December 31, 2020.
2. You don't have to purchase all boxes at once. Multiple invoices that reach the total qualifying boxes purchased should be submitted at one time using a single claim form.
3. All invoices must be received no later than January 31, 2021.
4. Packing slips are NOT acceptable.
5. Free offers are limited to 10 free boxes (*needles and syringes*) per customer location.
6. Each free goods redemption must meet the minimum required product purchase by catalogue number. Purchases are not limited to one catalogue number; however, the free item must have the same catalogue number as one of the qualifying purchased items.
7. The promotion is limited to catalogue numbers listed on the reverse side.
8. Wholesalers, distributors, acute care and customers accessing group purchasing contracts are not eligible for this program.

How to redeem your free product:

1. Fill out the information on this form.
2. Submit form and invoice(s) to CanadaFlu_Program@bd.com
3. BD will mail your free product. Customers will receive their product within 2–3 weeks.

Name: _____ Practice facility: _____

Address: _____

City: _____ Province: _____ Postal code: _____

Telephone: _____ Fax: _____ Email: _____

Distributor name: _____

☐ Yes, I would like to receive product information from BD by email to: _____@_____.

Protecting your privacy is very important to us. We will only use this email address to send you promotional material on our products. We will not disclose your personal information to any third parties. You will be able to opt out from this mailing list at any time.

Requested number of free boxes or cases:*

Number of free boxes or cases: _____ Product code: _____

*Free product catalogue number must be the same as the catalogue number for the purchased items.

BD–Canada, Mississauga, ON, L5N 0B3
905.288.6000

bd.com

BD, the BD Logo, Eclipse, Integra, Luer-Lok and SafetyGlide are trademarks of Becton, Dickinson and Company or its affiliates.
© 2020 BD. All rights reserved. (4928/0620)

